


 single use  
sterile packed



nxt™ single-use non-stick bipolar forceps



nxt™ non stick single-use bipolar forceps


## Superior non-stick technology - 100% made in Germany

- Solid silver tips for a superior non-stick effect (compared to silver coated tips)
- Slim forceps profile, thin tines and excellent visualization
- Can be used with standard bipolar cables
- Optimal balance and weight
- High stability (no scissoring)
- Precision coagulation
- True tip alignment
- Premium look & feel

*nxt™ Non-Stick Single-Use Bipolarpinzetten*

## *Überragende Non-Stick Eigenschaften - 100% made in Germany*

- *Solide Silberspitzen für einen überlegenen Non-Stick Effekt (im Vergleich zu beschichteten Pinzetten)*
- *Schmales Pinzettenprofil, dünne Schenkel, exzellente Sicht*
- *Verwendbar mit Standard Bipolarkabeln*
- *Optimale Balance und Gewicht*
- *Hohe Stabilität (kein Kreuzen der Pinzettenspitzen)*
- *Präzise Koagulation*
- *Exakt ausgerichtete Spitzen*
- *Exzellente Verarbeitung und Produktqualität*


nxt™  
Non-Stick Technology vs. **standard** bipolar forceps





wide handle

## single-use classic non-stick




**Ø 0.4 x 6 mm micro tip straight**

**standard series / working / total length 50 / 160 mm**

**BIP047USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™




**Ø 0.7 x 6 mm micro tip straight**

**standard series / working / total length 50 / 160 mm**

**BIP048USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™



**Ø 1.5 x 8 mm standard tip straight**

**standard series / working / total length 50 / 160 mm**


**BIP048-1.5USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™



wide handle

## single-use classic non-stick




Ø 0.4 x 6 mm micro tip straight

**standard series / working / total length 70 / 175 mm**

**BIP050USSU** us pin connector

alternative to Spetzler® Malis™ 6760-180-005 or VersaTru™ 9007050ST




Ø 0.7 x 6 mm micro tip straight

**standard series / working / total length 70 / 175 mm**

**BIP051USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™




Ø 1.0 x 8 mm standard tip straight

**standard series / working / total length 70 / 175 mm**

**BIP052USSU** us pin connector

alternative to Spetzler® Malis™ 6760-180-010 or VersaTru™ 9007100ST



Ø 1.5 x 8 mm standard tip straight

**standard series / working / total length 70 / 175 mm**


**BIP052-1.5USSU** us pin connector

alternative to Spetzler® Malis™ 6760-180-015 or VersaTru™ 9007150ST



wide handle

## single-use classic non-stick




### Ø 0.2 x 6 mm micro tip straight

**standard series / working / total length** 95 / 205 mm

**BIP053USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™




### Ø 0.4 x 6 mm micro tip straight

**standard series / working / total length** 95 / 205 mm

**BIP054USSU** us pin connector

alternative to Spetzler® Malis™ 6760(80)-200-005 or VersaTru™ 9008050ST(SL)




### Ø 1.0 x 8 mm standard tip straight

**standard series / working / total length** 95 / 205 mm

**BIP056USSU** us pin connector

alternative to Spetzler® Malis™ 6760(80)-200-010 or VersaTru™ 9008100ST(SL)




### Ø 1.2 x 8 mm standard tip straight

**standard series / working / total length** 95 / 205 mm

**BIP057USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™




### Ø 1.5 x 8 mm standard tip straight

**standard series / working / total length** 95 / 205 mm

**BIP056-1.5USSU** us pin connector

alternative to Spetzler® Malis™ 6760(80)-200-015 or VersaTru™ 9008150ST(SL)




### Ø 2.0 x 8 mm standard tip straight

**standard series / working / total length** 95 / 205 mm

**BIP058USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™




### Ø 1.2 x 8 mm standard tip 45° down

**standard series / working / total length** 95 / 205 mm

**BIP059USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™



### Ø 1.2 x 8 mm standard tip 45° up

**standard series / working / total length** 95 / 205 mm


**BIP060USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™



wide handle

## single-use classic non-stick




### Ø 0.4 x 6 mm micro tip straight

**standard series / working / total length** 130 / 235 mm

**BIP062USSU** us pin connector

alternative to Spetzler® Malis™ 6760(80)-230-005 or VersaTru™ 9009050ST(SL)




### Ø 1.0 x 8 mm standard tip straight

**standard series / working / total length** 130 / 235 mm

**BIP064USSU** us pin connector

alternative to Spetzler® Malis™ 6760(80)-230-010 or VersaTru™ 9009100ST(SL)




### Ø 1.2 x 8 mm standard tip straight

**standard series / working / total length** 130 / 235 mm

**BIP065USSU** us pin connector

alternative to Spetzler® Malis™ or VersaTru™



### Ø 1.5 x 8 mm standard tip straight

**standard series / working / total length** 130 / 235 mm


**BIP064-1.5USSU** us pin connector

alternative to Spetzler® Malis™ 6760(80)-230-015 or VersaTru™ 9009150ST(SL)




wide handle

## single-use irrigation non-stick



**Ø 0.5 x 6 mm micro tip straight**  
**standard series / working / total length 95 / 205 mm**  
**BIP002US-IRSU** us pin connector  
 alternative to Spetzler® Malis™ 6770-200-005



**Ø 1.0 x 8 mm standard tip straight**  
**standard series / working / total length 95 / 205 mm**  
**BIP003US-IRSU** us pin connector  
 alternative to Spetzler® Malis™ 6770-200-010







wide handle

1:1

## single-use irrigation non-stick



**Ø 0.5 x 6 mm micro tip straight**  
**standard series / working / total length 130 / 235 mm**  
**BIP004US-IRSU** us pin connector  
 alternative to Spetzler® Malis™ 6770-230-005



**Ø 1.0 x 8 mm standard tip straight**  
**standard series / working / total length 130 / 235 mm**  
**BIP005US-IRSU** us pin connector  
 alternative to Spetzler® Malis™ 6770-230-010





adeor® is a registered trademark of adeor medical AG.

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing device are reserved and require our written approval.

Copyright 2021 adeor medical AG

01-1121-D/E

Printed in Germany

November 2021

Spetzler and Malis are trademarks of Stryker Inc.

VersaTRU is a trademark of Johnson & Johnson Inc.

adeor is not affiliated with any of the above.



## HEADQUARTERS

adeor medical AG

Martinshof 5

D-83626 Valley

Germany

Phone +49 (0) 8024 47717 0

Fax +49 (0) 8024 47717 99

eMail [office@adeor.com](mailto:office@adeor.com)

web [www.adeor.com](http://www.adeor.com)

Printed on recycled paper.

Gedruckt auf Recyclingpapier.

