


External CSF Drainage Catheters

Sophysa has developed specific ventricular and lumbar catheters, in radiopaque silicone elastomer. They are meant to be used with external CSF drainage systems.



Ventricular catheters

Reference	Designation
External Ventricular Catheter (Medium Lumen) • Ventricular catheter in transparent silicone with a radiopaque stripe:	
	<p>DE-104</p> <ul style="list-style-type: none"> • 29 cm length, 2 mm ID, 3 mm OD, 1 mm diameter holes • Multiperforated over a distance of 23 mm from the proximal tip • Numerical markings from 3 to 28 cm from the proximal tip <p>Supplied with:</p> <ul style="list-style-type: none"> • A female Luer-Lock connector • Introducing stylet • Trocar for tunnelling • Suture tab
External Ventricular Catheter (Large Lumen) • Ventricular catheter in radiopaque silicone elastomer:	
	<p>FVPC4</p> <ul style="list-style-type: none"> • 35 cm length, 2.3 mm ID, 4 mm OD, 2 mm diameter holes • Multiperforated over a distance of 23 mm from the proximal tip • Depth markings at 5, 7.5, 10, 15 and 20 cm from the proximal tip <p>Supplied with:</p> <ul style="list-style-type: none"> • A female Luer-Lock connector • Introducing stylet • Trocar for tunnelling • Suture tab at 25 cm from the proximal end
External Ventricular Catheter • Ventricular catheter in transparent silicone with a radiopaque stripe:	
	<p>DE-106</p> <ul style="list-style-type: none"> • 29 cm length, 1.5 mm ID, 3.0 mm OD, 1 mm diameter holes • Multiperforated over a distance of 23 mm from the proximal tip • Depth markings at 5, 10, 15 cm from the proximal tip <p>Supplied with:</p> <ul style="list-style-type: none"> • Female Luer-Lock connector • Introducing stylet • Trocar for tunnelling • Suture tab


External CSF Drainage | Catheters



DE-109

Lumbar catheters

Product Code	Description
--------------	-------------



External Lumbar Catheter • Radiopaque silicone lumbar

DE-108

- 90 cm length, 0.76 mm ID, 1.60 mm OD, 0.5 mm diameter holes
- Multiperforated over a distance of 23 mm from the proximal tip
- Depth markings at 11, 16, 21, 26 cm from the proximal tip
- With an opened-end tip

Supplied with:

- Tuohy needle (14 G), 9 cm length
- Female Luer-Lock connector (compression hub)
- Suture tab

External Lumbar Catheter with Introducer Guide • Radiopaque silicone lumbar catheter:

DE-109

- 90 cm length, 0.76 mm ID, 1.60 mm OD, 0.5 mm diameter holes
- Multiperforated over a distance of 23 mm from the proximal tip
- Depth markings at 11, 16, 21, 26 cm from the proximal tip
- With a closed-end tip

Supplied with:

- Tuohy needle (14 G), 9 cm length
- Female Luer-Lock connector (compression hub)
- Suture tab
- Introducer guide



www.sophysa.us

Sophysa
 5, rue Guy Moquet
 91400 Orsay - France
 Tel.: +33 (0)1 69 35 35 00
 Fax: +33 (0)1 69 35 36 90
 contact@sophysa.com

Sophysa USA Inc.
 503 E Summit Street, Suite 5
 Crown Point, IN 46307 - USA
 Tel.: +1 (219) 663-7711
 Fax : +1 (219) 663-7741
 contact@sophysa.us

Distribution: